



CARESCAPE™ R860

Invasive Modes of Ventilation

March 19, 2020

Objectives

By the end of this course, you should be able to:

- Describe Compliance, Elastance and Airway resistance
- Categorize the CARESCAPE R860 ventilation modes
- Describe tube compensation, leak compensation and trigger compensation
- Identify and describe the CARESCAPE R860 modes of ventilation
- Discuss advantages and disadvantages for Volume Control, Pressure Control, Pressure Regulated Volume Control, and Airway Pressure Release Ventilation



Mechanics of Ventilation



Lung Compliance

Compliance = $\Delta\text{Volume}/\Delta\text{Pressure}$ (mL/cmH₂O)

Change in volume over the change in pressure

- A measure of the ease of expansion of the lungs and thorax, determined by pulmonary volume and elasticity.
- A high degree of compliance indicates a loss of elastic recoil of the lungs, as in old age or emphysema.
- Decreased compliance means that a greater change in pressure is needed for a given change in volume, as in atelectasis, edema, fibrosis, pneumonia, or absence of surfactant.

Static Compliance = Exhaled Tidal Volume (V_{te})/Plateau Pressure (P_{plat}) – Positive End Expiratory Pressure (PEEP)

Dynamic Compliance = Exhaled Tidal Volume (V_{te})/Peak Inspiratory Pressure (PIP) – Positive End Expiratory Pressure (PEEP)

Normal adult compliance 40-70 ml/cmH₂O

In children about 1ml/cmH₂O/kg



Elastance

$$\text{Elastance} = \Delta\text{Pressure} / \Delta\text{Volume (cmH}_2\text{O/mL)}$$


Change in pressure over the change in volume

- A measure of the tendency of something to recoil toward its original dimensions upon removal of a distending or compressing force.
- Compliance and elastance are inversely related.
 - If compliance increases then elastance decreases
 - If compliance decreases then elastance increases




Compliance and Elastance

Compliance



Alveolar Cavity



<https://www.visiblebody.com/blog/anatomy-and-physiology-gas-exchange>



Airway Resistance

$$R_{aw} = \Delta \text{ Pressure} / \text{Flow (cmH}_2\text{O/L/sec)}$$

Change in pressure over flow

- Airway resistance is the friction caused by the movement of air throughout the respiratory system

- Types of flow:
 - Laminar flow: smooth, even non-tumbling flow
 - Turbulent flow: rough, tumbling uneven flow pattern
 - The pressure gradient necessary to maintain turbulent flow is much higher than that necessary to maintain laminar flow.
 - Tracheobronchial flow: is a combination of laminar and turbulent flow which is maintained throughout the respiratory system

- Airway resistance decreases with increased airway diameter, bronchodilation, laminar flow and increase in lung volume

- Airway resistance increases with decreased airway diameter, bronchoconstriction, turbulent flow and decrease in lung volume

Normal airway resistance is 0.5-2.5cmH₂O/L/sec at a flow rate of 0.5 L/sec



Airway Resistance



Raw is normal or decreases



When the Raw increases



Modes of Ventilation Categories



Modes of Ventilation: Introduction and Overview

Define categories of ventilation

Control Modes:

Positive pressure ventilation in which the ventilator is in control mode, with its cycle entirely controlled by the apparatus and not influenced by the patient's efforts at spontaneous ventilation.

Synchronized Modes:

Synchronized Intermittent Mechanical Ventilation is a variation of IMV, in which the ventilator breaths are synchronized with patient inspiratory effort, with added pressure support.

Support Modes:

The patient initiates every breath and the ventilator delivers support with the preset pressure value. With support from the ventilator, the patient also regulates his own respiratory rate and tidal volume.



Modes of Ventilation

Ventilator mode can be defined as a set of operating characteristics that control how the ventilator functions.

Operating mode can be described by:

- The way a ventilator is triggered into inspiration and cycled into expiration.

- What variables are limited during inspiration.

- Whether or not the mode allows mandatory, spontaneous, or supported breaths.



Modes of Ventilation: Control Modes

Control Modes: Introduction and Overview

Each breath is initiated, limited and terminated by the ventilator.

Patients can breathe spontaneously between control breaths, but the ventilator does not respond to the spontaneous effort.

Mode	Primary Settings	Inspiratory Flow Pattern	Breath Timing	Patient Synchrony	Safety
AC/VC	Tidal Volume FiO2 PEEP Flow	Constant	Rate I:E, T _{insp} or T _{pause} Insp Pause	Insp Trigger Bias Flow	P _{limit} P _{max}
AC/PC	Inspiratory Pressure FiO2 PEEP	Decelerating	Rate I:E or T _{insp}	Insp Trigger Bias Flow Rise Time	P _{max}
AC/PRVC	Tidal Volume FiO2 PEEP	Decelerating	Rate I:E or T _{insp}	Insp Trigger Bias Flow Rise Time	P _{max} P _{min}



Modes of Ventilation: Synchronized Modes

Synchronized Modes: Introduction and Overview

Synchronized Intermittent Mechanical Ventilation in which the ventilator breaths are synchronized with patient inspiratory effort, with added pressure support.

Mode	Primary Settings	Inspiratory Flow Pattern	Breath Timing	Patient Synchrony	Safety
SIMV VC	Tidal Volume Flow FiO2 PEEP PS	Constant	Rate T _{insp} or T _{pause} Insp Pause	Insp Trigger Exp Trigger Bias Flow PS Rise Time	P _{limit} P _{max}
SIMV PC	Inspiratory Pressure FiO2 PEEP PS	Decelerating	Rate T _{insp}	Insp Trigger Exp Trigger Bias Flow Rise Time PS Rise Time	P _{max}
SIMV PRVC BiLevel VG	Tidal Volume FiO2 PEEP PS	Decelerating	Rate T _{insp}	Insp Trigger Exp Trigger Bias Flow Rise Time PS Rise Time	P _{max} P _{min}
BiLevel	Inspiratory Pressure FiO2 PEEP PS	Decelerating	Rate T _{insp}	Insp Trigger Exp Trigger Bias Flow Rise Time PS Rise Time	P _{max}



Modes of Ventilation: Support Modes

Support Modes: Introduction and Overview

The ventilator supplies pressure support in response to the spontaneous breathing with no set rate; pressure support can also be added to SIMV modes of ventilation.

The patient must be spontaneously breathing and the ventilator must recognize and respond to the spontaneous effort, based on the patient's inspiratory flow.

Mode	Primary Setting	Inspiratory Flow Pattern	Breath Timing	Patient Synchrony	Safety
CPAP/PS	PEEP PS FiO2	Decelerating	Patient Controlled	Insp Trigger Exp Trigger Bias Flow PS Rise Time	Pmax Minimum Rate Backup P _{insp} Backup T _{insp}
VS	Tidal Volume FiO2 PEEP	Decelerating	Patient Controlled	T _{supp} Insp Trigger Exp Trigger Bias Flow PS Rise Time	Pmax Pmin Minimum Rate Backup T _{insp}
APRV*	Phigh Plow FiO2	Decelerating	Thigh Tlow	Insp Trigger Bias Flow Rise Time	Pmax

* Indicates these modes could be either control or spontaneous modes of ventilation depending on patient effort



Ventilation Mode Features:

Tube Compensation

Leak Compensation

Trigger Compensation



Tube Compensation

- To set Tube compensation, a Tube Type and Tube Diameter must be set in the New Patient or Current Patient menu
 - The options for tube compensation are:
 - Endotrach
 - Trach
 - ---
 - When --- is selected, the ventilator will not compensate for tube resistance

- Provides additional pressure to compensate for the difference between the lung pressure and breathing circuit pressure during the inspiratory phase of pressure controlled and pressure-supported breaths
 - Can be used to offset all or a percentage of the additional resistive pressure created by the endotracheal tube

NOTE: Tube compensation increases the pressure delivered to the patient. The pressure delivered with tube compensation is limited to $P_{max} - 5 \text{ cmH}_2\text{O}$. Make sure that P_{max} is set appropriately for the patient when using tube compensation.



Leak Compensation

- When leak compensation is selected, a general message will show leak compensation is on
- When the ventilator detects a leak in the breathing circuit and leak compensation is active, the ventilator will respond in the following ways:
 - Flow and volume waveforms and measured volume data are adjusted to account for leaks
- The ventilator will adjust the tidal volume delivered to compensate for leaks in the following volume controlled modes:
 - A/C VC
 - A/C PRVC
 - SIMV VC
 - SIMV PRVC
 - BiLevel VG
 - VS
- The maximum tidal volume adjustment depends on the patient type:
 - Adult – 25% of the set tidal volume
 - Pediatric – 100% of the set tidal volume or 100ml, whichever is less
 - Neonatal – 100% of the set tidal volume

NOTE: The exhaled volume of the patient can differ from the measured exhaled volume due to leaks



Trigger Compensation

- Adjusts the flow trigger to compensate for leaks
 - Leaks can cause the ventilator to initiate breath automatically (auto-triggering)
- Trigger compensation reduces the need to manually adjust the inspiratory trigger setting to prevent auto-triggering.




Modes of Ventilation- Control Modes




Assist Control Volume Control (A/C VC)

- The ventilator delivers mechanical breaths of the set tidal volume



the set respiratory rate
pressure required to deliver the tidal volume
patient's lung compliance and resistance
able to synchronize mechanical breath to
patient's efforts and to allow triggering of
breaths.
patient can initiate spontaneous breaths at the
beginning of the expiratory phase
ventilator delivers an inspiratory flow based on the set
respiratory time and T_{pause}.
Flow is maintained during the inspiratory phase
until the pressure is below the pressure limit
If the pressure limit is reached, the gas flow is reduced to
the pressure limit level for the remainder of the
breath.
The ventilator delivers the set tidal volume and adjusts the
inspiratory flow as needed to maintain the set tidal
volume for patient breaths

Assist Control Volume Control (A/C VC)



1. Airway pressure (Paw) waveform
2. Inspiratory time (T_{insp})
3. Inspiratory pause (T_{pause})
4. Expiratory time (T_{exp})
5. PEEP
6. Flow waveform
7. Tidal volume (V_T)



Assist Control Pressure Control (A/C PC)

- The ventilator delivers mechanical breaths at the set inspiratory time at intervals based on the


...vered depends on the patient's

...le to synchronize mechanical breath to
...ous efforts and to allow triggering of
...breaths.


...nt can initiate spontaneous breaths at
...ring the expiratory phase
...urizes the circuit to the set

...atient decreases after the pressure
...ssure setting

...ses to maintain the set pressure for the
...y time



Assist Control Pressure Control (A/C PC)




1. Airway pressure (Paw) waveform
2. Inspiratory time (T_{insp})
3. Expiratory time (T_{exp})
4. Inspiratory pressure (P_{insp})
5. PEEP
6. Flow waveform
7. Tidal volume (V_T)




Assist Control Pressure Regulated Volume Control (A/C PRVC)

- The ventilator delivers mechanical breaths of the set tidal volume at intervals based on the set respiratory rate. For each breath, the ventilator adjusts the inspiratory pressure to use the pressure required to deliver the tidal volume. Ventilation settings may be different if breath timing have been changed. Due to the patient's lung compliance, the ventilator provides pressure-controlled ventilation for 10 seconds or 2 breaths, whichever is longer when the mode is initiated. Due to the patient's lung compliance, the inspiratory pressure is established for subsequent breaths. When adjusting the inspiratory pressure the following range is used:
 - PEEP + Pmin
 - Pmax - 2 cmH₂OThe difference in inspiratory pressure between breaths should not exceed +/- 3 cmH₂O. This feature is available to synchronize mechanical breaths with the patient's spontaneous efforts and to allow triggering of mechanical breaths. When the patient can initiate spontaneous breaths, the PEEP level during the expiratory phase



Assist Control Pressure Regulated Volume Control (A/C PRVC)



1. Airway pressure (Paw) waveform
2. Inspiratory time (T_{insp})
3. Expiratory time (T_{exp})
4. Variable pressure to deliver set TV
5. PEEP
6. Flow waveform
7. Tidal volume (VT)



Modes of Ventilation-

Synchronized Modes




Synchronized Intermittent Mandatory Ventilation Volume Control (SIMV VC)



provides synchronized mechanical breaths of a set tidal volume at intervals based on the set respiratory rate. Spontaneous efforts are delivered as spontaneous breaths. The pressure required to deliver the tidal volume is determined by the patient's lung compliance and resistance. The pressure limit settings may be different if breath timing (e.g., I-time and Flow) have been changed. The ventilator monitors delivered tidal volume and adjusts respiratory flow as needed to maintain the set tidal volume for subsequent breaths.

Synchronized Intermittent Mandatory Ventilation Volume Control (SIMV VC)



1. Airway pressure (Paw) waveform
2. Inspiratory time (T_{insp})
3. Inspiratory pause (T_{pause})
4. Spontaneous breathing period
5. Pressure-supported breath
6. Flow waveform
7. Tidal volume (V_T)
8. PEEP
9. Trigger window




Synchronized Intermittent Mandatory Ventilation Pressure Control (SIMV PC)



- The ventilator delivers synchronized mechanical breaths at a set pressure level for a set inspiratory time at the set respiratory rate. All other breaths are delivered as pressure-supported
- The time delivered depends on the patient's breathing effort
- No patient effort is available
- The ventilator pressurizes the circuit to the set pressure
- When the patient's effort decreases after the pressure is reached, the pressure setting remains constant
- When the patient's effort decreases to maintain the set pressure for the set inspiratory time


Synchronized Intermittent Mandatory Ventilation Pressure Control (SIMV PC)



1. Airway pressure (Paw) waveform
2. Inspiratory time (T_{insp})
3. Spontaneous breathing time
4. Pressure-supported breath
5. Inspiratory pressure (P_{insp})
6. Flow waveform
7. Trigger window
8. PEEP
9. Tidal Volume (V_T)




Synchronized Intermittent Mandatory Ventilation Pressure Regulated Volume Control (SIMV PRVC)



The ventilator delivers synchronized mechanical breaths of the intervals based on the set respiratory rate. At the start of each breath, the ventilator adjusts the pressure to use the lowest pressure required to achieve the set tidal volume. All other spontaneous efforts are pressure-supported breaths. The pressure support settings may be different if breath timing has been changed. If pressure support is available, the ventilator will provide controlled ventilation for 10 seconds or 2 breaths, whichever is longer when the mode is initiated. If the patient's lung compliance, the inspiratory pressure is established for subsequent breaths. The following settings are used:
- PEEP + Pmin
- PEEP max - 2 cmH2O
- The difference in inspiratory pressure between breaths does not exceed 7.5 cmH2O

Synchronized Intermittent Mandatory Ventilation Pressure Regulated Volume Control (SIMV PRVC)



1. Airway pressure (Paw) waveform
2. Inspiratory time (T_{insp})
3. Spontaneous breathing time
4. Variable pressure
5. PEEP
6. Flow waveform
7. Tidal volume (V_T)
8. Pressure supported breath
9. Trigger window




BiLevel Airway Pressure Ventilation Volume Guaranteed (BiLevel VG)

- The ventilator alternates between a set PEEP and the minimum pressure to deliver the set tidal volume based on the set rate and inspiratory time.



initiates a breath at the PEEP level, a pressure
 with at the PS settings is delivered.
 patient's lung compliance, the ventilator
 controlled ventilation for 10 seconds or 2
 whichever is longer when the mode is initiated.
 patient's lung compliance, the inspiratory
 established for subsequent breaths.
 the inspiratory pressure the following
 is used:
 $P + P_{min}$
 $P_{max} - 2 \text{ cmH}_2\text{O}$
 an inspiratory pressure between breath does
 $3 \text{ cmH}_2\text{O}$
 pressure alarm is active for the current
 breath's pressure target is $0.5 \text{ cmH}_2\text{O}$ lower
 is also available

BiLevel Airway Pressure Ventilation Volume Guaranteed (BiLevel VG)




1. Airway pressure (Paw) waveform
2. Inspiratory time (T_{insp})
3. Spontaneous breathing time
4. Variable pressure
5. PEEP
6. Flow waveform
7. Tidal volume (VT)
8. Pressure supported breath
9. Trigger window




BiLevel Airway Pressure Ventilation (BiLevel)*

- The ventilator alternates between the set PEEP level and the pressure level based on the set rate and

breath spontaneously at either level. If a breath is initiated at the PEEP level, a pressure-supported breath is initiated during the high pressure level of inspiratory pressure provided by P_{insp} settings. If a breath is initiated near the end of T_{high}, the ventilator provides the pressure support to support the breath. If a breath is initiated near the end of T_{high}, the ventilator provides no pressure support. If a breath is initiated near the end of T_{high}, the ventilator provides the pressure support to support the breath. If a breath is initiated near the end of T_{high}, the ventilator provides no pressure support. If a breath is initiated near the end of T_{high}, the ventilator provides the pressure support to support the breath. If a breath is initiated near the end of T_{high}, the ventilator provides no pressure support.



BiLevel Airway Pressure Ventilation (BiLevel)*



1. Airway pressure (Paw) waveform
2. Tinsp
3. Exp time
4. Pressure Support (PS)
5. Pinsp
6. PEEP
7. Flow waveform
8. Tidal Volume (VT)



Modes of Ventilation- Support Modes



Continuous Positive Airway Pressure /Pressure Support (CPAP/PS)

- The ventilator maintains a PEEP level and provides pressure support.

used on spontaneously breathing patients who take spontaneous breaths and determines the timing, and tidal volume.


When Rate is set, the ventilator will deliver a mechanical breath if the patient's respiratory rate is less than the minimum rate.

A backup breath will be delivered at the Backup Pressure setting for the time duration of the Backup Time setting.

Pressure Support is also available



Continuous Positive Airway Pressure /Pressure Support (CPAP/PS)




1. Airway pressure (Paw) waveform
2. Pressure support (PS)
3. Inspiratory time (Backup T_{insp})
4. PEEP
5. Flow waveform
6. Backup P_{insp}
7. Minimum rate backup breath
8. Tidal volume (TV)



Volume Support (VS)

- The patient initiates spontaneous breaths and determines respiratory rate and timing. The ventilator maintains a PEEP level and provides support to deliver the set tidal volume.
 - Intended for spontaneously breathing patients



The ventilator adjusts the inspiratory pressure to maintain the pressure required to deliver the tidal volume. If the patient's lung compliance is low, the ventilator delivers a backup ventilation for 10 seconds or 2 breath periods, whichever is longer, when the mode is initiated.

Based on the patient's lung compliance, the inspiratory pressure is adjusted for subsequent breaths.

The inspiratory pressure the following

is set:

P_{min}

2 cmH₂O

The inspiratory pressure between breaths does

not

change. If the rate is set, the ventilator will deliver a

mechanical breath if the patient's

respiratory rate is less than the minimum rate.

A backup breath will be delivered at the Backup PRVC


if the duration of the Backup Inspiratory Time

expires and a pressure alarm is active for the current breath,

the pressure target is 0.5 cmH₂O lower

Backup ventilation is also available

Volume Support (VS)



1. Airway pressure (Paw) waveform
2. Spontaneous inspiratory time
3. Spontaneous breathing period
4. Variable pressure
5. PEEP
6. Flow waveform
7. Tidal Volume (VT)
8. Inspiratory time (Backup Tinsp)
9. Backup P_{Insp}
10. Minimum rate backup breath




Airway Pressure Release Ventilation (APRV)*

- The ventilator alternates between a set high (P_{high}) and a set low (P_{low}) pressure level. On spontaneously breathing patients, the ventilator delivers the set (P_{high}) pressure for the set time. The ventilator then delivers the set (P_{low}) pressure for the set time. Spontaneous breaths are allowed to occur at either pressure level.



Airway Pressure Release Ventilation (APRV)*



1. Airway pressure (Paw) waveform
2. Thigh
3. Tlow
4. Phigh
5. Plow
6. Flow waveform



Modes of Ventilation Advantages and Disadvantages



Advantages and Disadvantages of Volume Control and Pressure Control

	Advantages	Disadvantages
Volume Control	<ul style="list-style-type: none">• Constant tidal volume• Consistent alveolar ventilation• Easily identify changes in PIP and Pplat as respiratory mechanics change	<ul style="list-style-type: none">• Constant flow rate• Increase in potential asynchronies• Varying pressures
Pressure Control	<ul style="list-style-type: none">• PIP and peak alveolar pressures are constant• Flow varies with patient demand	<ul style="list-style-type: none">• Varying tidal volumes



Advantages and Disadvantages of Pressure Regulated Volume Control

	Advantages	Disadvantages
Pressure Regulated Volume Control	<ul style="list-style-type: none">• Targeted tidal volume• Pressure automatically adjusts based on lung compliance and airway resistance• Decelerating waveform• Variable inspiratory flow to meet patient's demand	<ul style="list-style-type: none">• Pressure adjusts based on the tidal volume of the last breath• Asynchronies may occur with variable patient effort



Advantages and Disadvantages of Airway Pressure Release Ventilation

	Advantages	Disadvantages
Airway Pressure Release Ventilation	<ul style="list-style-type: none">• Uses “Open lung” concept<ul style="list-style-type: none">• Maximize and maintain alveolar recruitment• Improve Oxygenation• Potential lung protective effect• Preservation of spontaneous breathing<ul style="list-style-type: none">• Less need for sedation and neuromuscular blocking agents• Better ventilation to dependent lung regions• Better cardiac filling with spontaneous breathing• Reduce the risk of ventilator induced diaphragmatic dysfunction	<ul style="list-style-type: none">• Increased work of breathing and oxygen consumption with spontaneous breathing• May create asynchrony and discomfort• Potential risks for volutrauma<ul style="list-style-type: none">• Large tidal volume swings with spontaneous effort• Increased transpulmonary pressures• Greater need for clinical trials to demonstrate better clinical outcomes over conventional ventilation



Daoud EG, Farag HL, Chatburn RL. Airway Pressure Release Ventilation: What Do We Know? *Respiratory Care*. 2012;57(2):282. doi:[10.4187/respcare.01238](https://doi.org/10.4187/respcare.01238)

Mireles-Cabodevila E, Kacmarek RM. Should Airway Pressure Release Ventilation Be the Primary Mode in ARDS? *Respiratory Care*. 2016;61(6):761. doi:[10.4187/respcare.04653](https://doi.org/10.4187/respcare.04653)

Myers TR, Macintyre NR. Does Airway Pressure Release Ventilation Offer Important New Advantages in Mechanical Ventilator Support? *Respiratory Care*. 2007;52(4):452.

Conclusion

This concludes the CARESCAPE R860 Modes of Ventilation.

In this course, you learned about:

- Compliance, Elastance and Airway resistance
- CARESCAPE R860 ventilation modes categories
- Tube compensation, leak compensation and trigger compensation
- The CARESCAPE R860 modes of ventilation
- The advantages and disadvantages of volume control, pressure control, pressure regulated volume control and airway pressure release ventilation

Disclaimers

Always refer to device manufacturers user reference manual for specific application of your CARESCAPE R860. Use this information as guidance and each patient may require clinical decisions not covered in this information. Ensure proper clinically appropriate alarm limits are set and monitored.



