

# IHE Integration Statement

Vendor	Product Name	Version	Date
Image Diagnost International GmbH	MammoWorkstation	4.5	Mai 11, 2009

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profiles Implemented	Actors Implemented	Option Implemented
Mammography Image	Image Display	None
	Evidence Creator	None
	Print Composer	None
Scheduled Workflow (SWF)	Image Display	None
	Evidence Creator	None
Consistent Presentation of Images (CPI)	Image Display	None
	Evidence Creator	None
	Print Composer	None
Access to Radiology Information (ARI)	Image Display	None
	Report Reader	None
Consistent Time	Time Client	None
Evidence Documents	Image Display	None
	Evidence Creator	None
Simple Image and Numeric Report	Report Creator	None
	Report Reader	None

**Internet address for vendor's IHE information:**

<http://www.imagediagnost.de/conformance>

**Links to Standards Conformance Statements for the implementation :**

<http://www.imagediagnost.de/conformance>

**Links to general information on IHE:**

<http://www.ihe-europe.org>

**Image Diagnost international GmbH**

Balanstr. 57  
81541 München

Tel +49 (89) 450 366 6  
Fax +49 (89) 450 366 89  
[www.imagediagnost.de](http://www.imagediagnost.de)

