



Reusable Anesthesia Breathing Masks

Reusable Anesthesia Breathing Masks from GE Healthcare are tested and verified to function with GE Healthcare anesthesia machines. Anesthesia breathing masks are designed to channel gas mixtures to the patient for non-invasive ventilation via the mouth and nose. The products perform optimally with GE Healthcare anesthesia machine breathing systems when applicable.

FEATURES

- Clear face piece allows visual observation of patient status
- Compatible with all standard breathing systems with 22mm and 15mm connectors where applicable
- Available in sizes from neonatal through large adults
- Reusable masks are autoclavable up to 50 times

**Our equipment.
Our accessories.
Simply compatible.**

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- Simplified ordering process
- Fast delivery
- One contact point for expert advice
- Quality support at every step

REUSABLE ANESTHESIA MASK OFFERINGS

Silicone reusable anesthesia masks are designed to channel anesthesia gases or fresh gas mixtures to the patient for inhalation via the mouth and nose. Anesthesia masks are to be used in conjunction with anesthesia and ventilator breathing systems. They are reusable and fully autoclavable up to 50 cycles. The molded texture of the rim ensures maximum patient comfort and a tight seal, while the transparent face piece allows patient monitoring at all times.

PRODUCT SPECIFICATIONS

Material composition

Component	Material
Mask	Silicone
Connector	Silicone
Silicone Reusable Face Masks	Do not contain natural rubber latex.
Silicone Reusable Face Masks	Do not contain natural phthalate DEHP.

Product use

Silicone Reusable Anesthesia Face Masks are reusable devices. They can be reused up to 50 times.

Sterility

Silicone Reusable Anesthesia Face Masks are supplied non-sterile.

Storage

Store in a cool, dry place out of direct sunlight

Shelf life

Shelf life of 5 years from the date of manufacture. This is based on the stability of the devices' components and raw materials sourced. Expiry date is clearly marked on individual product pouch.

PART LIST

Part Number	Description	Quantity
2106570-011	Mask, reusable, size 0, 15mm connection	1
2106570-019	Mask, reusable, size 1, 15mm connection	1
2106570-020	Mask, reusable, size 2, 22mm connection	1
2106570-021	Mask, reusable, size 3, 22mm connection	1
2106570-022	Mask, reusable, size 4, 22mm connection	1
2106570-023	Mask, reusable, size 5, 22mm connection	1
2106570-029	Mask, reusable, size 0, 15mm connection	10
2106570-030	Mask, reusable, size 1, 15mm connection	10
2106570-031	Mask, reusable, size 2, 22mm connection	10
2106570-032	Mask, reusable, size 3, 22mm connection	10
2106570-033	Mask, reusable, size 4, 22mm connection	10
2106570-034	Mask, reusable, size 5, 22mm connection	10

Disposal considerations

Dispose as clinical waste, in accordance with hospital policy, local guidelines and regulations.

Packaging materials

Primary	Polybag – Polyethylene (PE)
Secondary	Carton / Box – Cardboard

Product may not be available in all countries and regions.

Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com

Distributed by:
Datex-Ohmeda, Inc.



Manufactured by:
Flexicare Medical Limited
Cynon Valley Business Park,
Mountain Ash, CF45 4ER, UK