



Intrauterine Pressure Catheters and Cables

Quality accessories. Easy to use.

In labor and delivery, patient comfort and safety are top priorities. Our Corometrics™ monitor family offers features you know and trust, including the latest advancements in maternal and fetal monitoring and high-quality monitoring accessories.

GE intrauterine pressure catheters and cables are designed to help you deliver high-quality care by providing accurate measurements.

Catheters

- Easy catheter insertion with a rounded transducer tip
- Flat catheter body to deflect away from tissue, minimizing trauma during insertion
- Fast liquid infusion with a flexible, clear amnioinfusion line with tethered cap, delivering gravity feed at 55 ml/mm
- Inaccurate readings minimized via splash-resistant connector

Cables

- Easy identification of cable issues with built-in cable check feature
- Convenient zero/re-zero with a single button and indicator light for re-zero confirmation

Our equipment. Our accessories. Simply compatible.

You can trust GE Healthcare accessories and consumables to keep your clinical devices working efficiently.

- One contact point for expert advice
- Simplified ordering process
- Fast delivery
- Quality support at every step

FEATURES

Product compatibility

Cables are compatible and verified with the following monitoring devices:

- Corometrics 174 fetal monitor
- Corometrics 259cx fetal monitor
- GE Mini Telemetry System

PRODUCT SPECIFICATIONS

Materials

Not made with natural rubber latex or DEHP
ROHS (Restriction of hazardous substances) compliant

Operating conditions

Recommended ambient temperature 10° to 35°C (50° to 95°F)

Storage conditions

Recommended temperature -10° to 50°C (14° to 122°F)

Cleaning solutions (cable only)

Ethanol	70-96% by volume
Isopropyl alcohol maximum	60% by weight
Glutaraldehyde	2%
Phenol	2%

The catheter is single use, and the cable is re-usable.

Compliance

Sold only in US and Canada

Product may not be available in all countries and regions.

Contact a GE Healthcare Representative for more information. Please visit www.gehealthcare.com/promotional-locations

PART LIST

Part Number	Description	Quantity
2108982-001	56300 Intrauterine pressure catheter, single use	10
2108998-001	56321 Intrauterine pressure catheter cable	1



Manufactured by:
Cardinal Health 200, LLC
7000 Cardinal Place
Dublin, OH 43017
USA